

Molecular virology

- Laboratory exercises -

VIRAL NUCLEIC ACIDS (NA) ISOLATION


Sandra Kolundžija, Ph.D.

Asst. Prof. Silvija Černi, Ph. D.

Lucija Nuskern Karaica, Ph. D.


Prof. Dijana Škorić, Ph. D.

VIRUSES – intracellular pathogens




Isolation of viral RNAs from plant tissue:

- Total RNAs (plant RNA + viral RNA; fraction of viral RNAs is small relative to plant RNAs)
 - viral dsRNAs (replicative form)
- **Isolation of genomic viral RNAs from previously partially purified viral particles** (increasing the fraction of viral RNAs in sample, i.e. separating it from plant RNAs)




VIRAL NUCLEIC ACIDS ISOLATION (RNA):



Tobacco mosaic virus (TMV)

ssRNA(+) genome, 6.3 - 6.5 kb
anisometric virus
5 % NA in the virus particle



Turnip yellow mosaic virus (TYMV)


ssRNA(+) genome, 6.3 kb
isometric virus
28 % NA in the virus particle

PURIFICATION OF VIRUS PARTICLES



HOMOGENISATION OF
INFECTED TISSUE
(+ ribonuclease inhibition)

ORGANIC SOLVENT
ADDITION




water phase (virus particles, leftover plant proteins, leftover ribosomes, low-molecular components)

clot/interphase (rough cellular debris, denatured proteins and nucleoproteins on surface)

organic solvent phase (pigments, membranous structures)

CLARIFICATION


ULTRACENTRIFUGATION
or alternative PHYSICO
CHEMICAL procedure
(addition of ammonium
sulphate + low speed
centrifugation)



partially purified virus solution

plant proteins

From here, we isolate viral NA.



TMV

E=3

$A_{260}/A_{280} = 1.19$


TYMV

E=9,6

$A_{260}/A_{280} = 1.51$


Traditional Phenol Extraction

Phenol, pH 4



RNA

Phenol, pH 8



DNA

Disadvantages:
Phenol – toxic, carcinogenic




→ Alternative methods

PERCHLORATE METHOD OF NA ISOLATION – specific steps

1. SDS (Sodium dodecyl sulphate) – anionic detergent

- amphiphilic molecule – it disrupts interactions between non-polar amino acids (AA) in proteins – protein denaturation (enhanced by heating)

$\text{NaC}_{12}\text{H}_{25}\text{SO}_4$ = Sodium dodecyl sulfate




Globular proteins + detergent molecules = micellae (good solubility, clear solution)

PERCHLORATE METHOD OF NA ISOLATION – specific steps

2. Sodium perchlorate (NaClO_4) – hygroscopic

- It extracts water from hydration shell of micellae – turbidity („flocks”) →
- Proteins separate from aqueous phase containing NA
- „Flocks” are lighter than the aqueous phase – they float



3. Separation of aqueous phase with NA – centrifuge

4. NA precipitation with **EtOH** (3 volumes for RNA, 2 volumes for DNA)

PERCHLORATE METHOD OF VIRAL NUCLEIC ACIDS (RNA) ISOLATION – protocol

+ 0.8 mL partially purified virus
+ 0.2 mL **SDS** (25%)

heat for 3 min at 60°C

denaturation of capsid proteins

+ 3 mL **NaClO₄** (8 M)

vortex for 1 min

separation of proteins and NA (turbidity)


- centrifuge for 10 min at 6000 rpm (4 °C)
- separate aqueous phase containing RNA (at the bottom of the tube), add 3V of cold **EtOH** and incubate for 20 min at -20 °C
- centrifuge for 20 min at 11000 rpm (4 °C)
- dry the precipitate and dissolve it in 25 µL **TE-buffer** (pH 7.5)

TE buffer

- 10 mM Tris + 1mM EDTA
- pH 7.5 for RNA, pH 8 for DNA
- nuclease inactivation – protects NA from degradation

PARTIAL DENATURATION OF NA (preparation for electrophoresis)

•Formamide



- Disrupts hydrogen bonds between base pairs → stabilisation of ssRNAs in electrophoretic gel
- Affects T_m (higher order structures dissociate at lower temp.)
- NA in semi-denatured state: better electrophoretic properties, stability, retains infectivity (possibility of downstream biological experiments)

PARTIAL DENATURATION OF NA – protocol

10 μ L nucleic acid isolate
+ 40 μ L **formamide**

heat for 2 min at 85 °C

add 10 μ L **coloured electrophoresis sample buffer (SB)**

ELECTROPHORESIS OF NUCLEIC ACIDS – protocol

Agarose gel preparation:

0.5% agarose	1% agarose	1.2% agarose
150 mg agarose	300 mg agarose	360 mg agarose
30 mL 1xTBE	30 mL 1x TBE	30 mL 1x TBE

Agarose = agar purified from agaropectin (charged groups disrupt electrophoresis – electro-osmosis)

Electrophoresis:

90 V, 40 min



Xylene cyanol FF
TAE: 4160bp
TBE: 3030bp

Bromophenol blue
TAE: 370bp
TBE: 220bp

Electrophoresis of tracking dye
in 6xGel Loading Dye, SDS +

Composition and role of SB (sample buffer, gel loading buffer, loading „dye”)

- dyes (xylene cyanol, bromphenol blue) → sample visibility (loading, run tracking)
- sample density increase (sucrose, glycerol, ficoll)
- water

ADVANTAGES AND DISADVANTAGES OF THE PERCHLORATE METHOD

+

- **Fast**
- **Non-toxic**
- **Excess of perchlorate in EtOH inhibits the precipitation of proteins leftover in the aqueous phase**

-

- **Samples are never as pure as after phenol purification**

ALTERNATE METHODS FOR NA ISOLATION

- **Phenol method – „classical method”;**
 - phenol – toxic, mutagen, harmful for environment;
 - method: slower, „finer” results
- **“Flash heating”**
 - opening of virus particle by heating (physical method)
 - several cycles of heating for 1 min at 95°C
 - risk of NA denaturation

SILVER STAINING

- Very sensitive method, developed for PAGE (Ag^+ does not bind strongly to polyacrylamide gels, but it does bind to agarose)

1. Fixation:


- Formaldehyde (reducing agent) – increase of contrast due to blocking of charged groups of agarose
- Simultaneous NA fixation (to enable diffusion) and complete denaturation

2. Washing:

- Rinsing of buffer, ions, denaturing agent
- Formaldehyde should be COMPLETELY rinsed off! (to avoid unspecific Ag^+ reduction in staining)

3. Silver impregnation:

- AgNO_3 – binds to gel and NA
(rinsing removes AgNO_3 unspecifically bound to gel)



Impregnation
with silver ion

Reduction of silver
ion to metallic silver

4. Developing:

- Developing solution – Ag^+ reduction to Ag^0 forms a precipitate (formaldehyde = reducing agent, sodium carbonate increases pH)

SILVER STAINING – protocol

Solution A (FIXATION + REDUCTION)

110 mL water
+ 14 mL formaldehyde

incubate gel for 20 min; 3 x 10 min rinsing with water

Solution B (SILVER IMPREGNATION)

250 mg AgNO₃
+ 125 mL water

incubate gel for 10 min in dark, 3 x short rinse with water (15")

Solution C (DEVELOPING SOLUTION)

6,867 g Na₂CO₃ x 10 H₂O
70 µL formaldehyde
+ 125 mL water

incubate gel until bands become visible

Solution D (FIXATIVE)


0,5 mL HNO₃
+ 125 mL water

incubate gel for several min

Store the gels in water in dark!



NA STAINING WITH NON-FLUORESCENT DYES

- Methylene blue
- Toluidine blue



NA STAINING WITH FLUORESCENT DYES


- Ethidium bromide (EtBr)
 - max A (EtBr) = 360 nm
 - max A (NA) = 260 nm – nucleotides excited by UV-light transfer energy to EtBr
- Acridine orange – differentiates DNA (green) and RNA (red)
- Commercially available fluorescent dyes (Sybr Green, Sybr Safe, GelRed, Gel Green, Stain G...)



MARKERS


1 kb DNA Ladder, 250 - 10,000 bp

0.25 – 10.0 kb (ng/0.5 µg)



1.0 % Agarose

**ΦX174 DNA/BsuRI (HaeIII)
Marker, 9, ready-to-use**



0.5 µg/lane, 8 cm length gel,
1X TBE, 5 V/cm, 1.5 h

RESULTS – GROUP 1

Electrophoresis:

1.2% gel, 1xTBE

40 min, 100 V

a) Silver staining


b) Ethidium bromide staining

a)

M TMV TYMV



0.25 – 10.0 kb (ng/0.5 µg)



1.0 % Agarose


Total RNA concentration after NA isolation:

	Group 1
TMV	390,1 µg/mL
TYMV	1273,6 µg/mL

Load : TMV ~ 1.5 - 2 µg per well
TYMV ~ 5.2 - 6.5 µg per well

b)

M TMV TYMV



viral genomic RNA

decomposed RNA

RESULTS – GROUP 2


Total RNA concentration after NA isolation:

	Group 2
TMV	241,4 µg/mL
TYMV	217,8 µg/mL


Load:

- 1 – TMV - 0.5 µg
- 2 – TMV - 0.75 µg
- 3 – TMV - 1 µg
- 4 – TYMV - 0.5 µg
- 5 – TYMV - 0.75 µg
- 6 – TYMV - 1 µg

a)
1% gel, 1xTBE, 40 min, 100 V
Silver staining



b)
0.5% gel, 1xTBE + 1 µL Serva DNA Stain G
~1 h, 100 V



RESULTS – GROUP 3

Total RNA concentration after NA isolation:

	Group 3
TMV	426,6 µg/mL
TYMV	83,1 µg/mL

Load:

- 1 – TMV - 1 µg
- 2 – TYMV - ?
- 3 – ?

a)

1% gel, 1xTBE, 40 min, 110 V
Silver staining

M 1 2 3 M9



b)

0.5% gel, 1xTBE + 1µL Serva DNA Stain G
~50 min, 110 V

M 1 2 3 M9

